



---

# Product Catalog

October 2021

---

[www.reachsurgical.com](http://www.reachsurgical.com)



## **EXPAND THE SCOPE OF SURGERY**

Reach Surgical, a Genesis Medtech company, is committed to serving high-quality products to clinicians and patients worldwide.

We deliver advanced technology and cost-effective minimally invasive surgical products to medical institutions and provide patients with reliable and economical treatment solutions.

## CONTENTS

Endoscopic Stapling	4-7
Open Stapling	8-9
Energy	10-11

## ENDO III Stapler



Series	Product Code	Handle Length (mm)	Max number of reloads	QTY/BOX
ENDO III 45mm Articulating Stapler	REC45A S	282	12	3
	REC45A	332		
	REC45A L	432		
ENDO III 60mm Articulating Stapler	REC60A S	298		
	REC60A	348		
	REC60A L	448		




## ENDO III Cartridge

Product Code	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Tissue Thickness	Number of staple rows	QTY/BOX
REC 45 GRA	Gray	2.0	0.75	Extra Thin	6	12
REC 60 GRA						
REC 45 WHT	White	2.5	1.0	Vascular		
REC 60 WHT						
REC 45 BLU	Blue	3.5	1.5	Regular		
REC 60 BLU						
REC 45 GLD	Gold	3.8	1.75	Regular/Thick		
REC 60 GLD						
REC 45 GRN	Green	4.1	2.0	Thick		
REC 60 GRN						
REC 45 BLK	Black	4.4	2.2	Extra Thick		
REC 60 BLK						

## iReach Magnum Powered Stapler

Product Code	Staple Line Length (mm)	Shaft Length (mm)	Instrument Length (mm)	Trocar Compatibility	Max number of reloads	QTY/BOX
IM60AS	60	298	546	12 mm	12	3
IM45AS	45	282	530			
IM60AM	60	348	596			
IM45AM	45	332	580			
IM60AL	60	448	696			
IM45AL	45	432	680			





## ENDO SRC Universal Stapler

Product Code	Shaft Length (mm)	Handle Length (mm)	Max number of reloads	QTY/BOX
ENDO SRC S	65	696	25	3
ENDO SRC	155	680		
ENDO SRC L	255	596		

## ENDO SRC Single Use Loading Unit



Product Code *	Staple Line Length (mm)	Cartridge Color	Convex Points	Open Staple Height (mm)	Closed Staple Height (mm)	Tissue Thickness	Trocar Compatibility	Number of staple rows	QTY/BOX
ENDO SRC4525 L ENDO SRC4525 R / ENDO SRC4525 BR	45	White	N/A	2.5	1.0	Thin	12 mm	6	6
ENDO SRC6025 L ENDO SRC6025 R / ENDO SRC6025 BR	60								
ENDO SRC4535 L ENDO SRC4535 R / ENDO SRC4535 BR	45	Blue	Yes	3.5	1.5	Regular			
ENDO SRC6035 L ENDO SRC6035 R / ENDO SRC6035 BR	60								
ENDO SRC6040 L ENDO SRC6040 R / ENDO SRC6040 BR	60	Gold	Yes	4.0	1.75	Regular/ Thick			
ENDO SRC4548 L ENDO SRC4548 R / ENDO SRC4548 BR	45	Green	Yes	4.8	2.0	Thick			
ENDO SRC6048 R	60								

\* L for Straight, R for Articulating, BR for Articulating with curved tip.

## ENDO AFT Single Use Loading Unit

Product Code *	Staple Line Length (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Tissue Thickness	Trocar Compatibility	Number of staple rows	QTY/BOX
ENDO AFT 45TNR ENDO AFT 45TNBR	45	Tan	2/2.5/3	0.75/1/1.25	Thin	12 mm	6	6
ENDO AFT 60TNR ENDO AFT 60TNBR	60							
ENDO AFT 45PLR ENDO AFT 45PLBR	45	Purple	3/3.5/4	1.25/1.5/1.75	Regular/ Thick			
ENDO AFT 60PLR ENDO AFT 60PLBR	60							
ENDO AFT 45BKR	45	Black	4/4.5/5	1.75/2/2.2	Thick	15 mm		
ENDO AFT 60BKR	60							

\* R for Articulating, BR for Articulating with curved tip.

## iReach SRC Powered Stapler

Product Code	Shaft Length (mm)	Handle Length (mm)	Max number of reloads	QTY/BOX
ENDO SRC PSZ/PS	72	326	25	3
ENDO SRC PZ/P	162	416		
ENDO SRC PLZ/PL	262	516		



## ACS-D Circular Stapler



Product Code	Color	Shaft Length	Outer Diameter (mm)	Lumen (mm)	Open Staple Height (mm)	Closed Staple Height (mm)	Number of Staples	Rows of Staples	QTY/BOX
ACS 21D	Gray	27 cm	21.5	12.5	3.5	1.5-1.8	16	2	3
ACS 25D	White		25.5	16	4.8	1.5-2.2	22		
ACS 28D	Blue		28.5	18.8			26		
ACS 31D	Green		31.5	22.7			30		
ACS 34D	Gold		34.5	25			32		

This product is not available for Germany

## PPH 33D

### Procedure Set for Hemorrhoid and Prolapse



Product Code	Cutting Area Diameter (mm)	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staples	Number of staple rows	QTY/BOX
PPH 33D	24.4	3.5	0.8 ~ 1.5	32	2	3

## RLS Linear Stapler



Product Code	Stapler Line (mm)	Preloaded with unit	Max Firings	QTY/BOX
RLS3035S	30	RLS3035L	8	3
RLS4535S	45	RLS4535L		
RLS6035S	60	RLS6035L		
RLS3048S	30	RLS3048L		
RLS4548S	45	RLS4548L		
RLS6048S	60	RLS6048L		

## RLS Linear Stapler Cartridge

Product Code	Stapler Line (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staples	Number of staple rows	QTY/BOX
RLS3035L	30	Blue	3.5	1.5	11	2	6
RLS4535L	45				15		
RLS6035L	60				21		
RLS3048L	30	Green	4.8	2.0	11		
RLS4548L	45				15		
RLS6048L	60				21		



## FLC Linear Cutter



Product Code	Stapler Line Length (mm)	Preloaded with cartridge	Max Firings	QTY/BOX
FLC6038S	60	FLC6038L	8	3
FLC6048S	60	FLC6048L		
FLC8038S	80	FLC8038L		
FLC8048S	80	FLC8048L		

## FLC Linear Cutter Cartridge




Product Code	Stapler Line Length (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staple rows	QTY/BOX
FLC6038L	60	Blue	3.8	1.5	4	6
FLC6048L	60	Green	4.8	2.0		
FLC8038L	80	Blue	3.8	1.5		
FLC8048L	80	Green	4.8	2.0		

## RLC Linear Cutter



Product Code	Stapler Line Length (mm)	Preloaded with cartridge	Max Firings	QTY/BOX
RLC6038S	60	RLC6038L	8	3
RLC6048S	60	RLC6048L		
RLC8038S	80	RLC8038L		
RLC8048S	80	RLC8048L		

## RLC Linear Cutter Cartridge



Product Code	Stapler Line Length (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staple rows	QTY/BOX
RLC6038L	60	Blue	3.8	1.5	4	6
RLC6048L	60	Green	4.8	2.0		
RLC8038L	80	Blue	3.8	1.5		
RLC8048L	80	Green	4.8	2.0		

## Sound Reach energy platform

Product Code	Description	QTY/BOX
CSUS8000	Generator	1
TRA6	Transducer	1
FSW2	Foot switch	1



## Sound Reach Shears

Product Code	Shaft Length (cm)	QTY/ BOX
CH14PD	14	3
CH23PD	23	3
CH36PD	36	3
CH45PD	45	3

This product is not available for Germany



## Contacts

**Global Headquarter**  
120 Xinxing Road  
West Zone, TEDA  
300462 Tianjin  
China

**T** +86-22-25326259  
**F** +86-22-66200915

## Importer



## Distributor

