


AndraBalloon® & Andra2Balloon
PTA-Catheter

Optimized for stent dilatation
Balloon diameter 10/14 – 32 mm
Balloon length 20 – 60 mm

AndraBalloon® & Andra2Balloon PTA-Catheter



- Optimized folding for better rewrapping and easy repassage
- Coaxial shaft design for fast inflation and deflation
- Optimized PTA-balloon for stent dilatation
- 2 Pt/Ir X-ray markers on the inner shaft

Order Codes AndraBalloon® PTA-Catheter

Order Code	Diameter/Length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
ABL 10/20-1200 ABL 10/30-1200 ABL 10/40-1200	10/20 10/30 10/40	400 (4)	1 000 (10)	7
ABL 12/20-1200 ABL 12/30-1200 ABL 12/40-1200	12/20 12/30 12/40	400 (4)	1 000 (10)	7
ABXL 14/30-1200 ABXL 14/40-1200 ABXL 14/60-1200	14/30 14/40 14/60	400 (4)	1 000 (10)	8
ABXL 16/30-1200 ABXL 16/40-1200 ABXL 16/60-1200	16/30 16/40 16/60	400 (4)	800 (8)	10
ABXL 18/30-1200 ABXL 18/40-1200 ABXL 18/60-1200	18/30 18/40 18/60	400 (4)	700 (7)	10

Packaging Unit: 1

Order Code	Diameter/Length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
ABXL 20/30-1 200 ABXL 20/40-1 200 ABXL 20/60-1 200	20/30 20/40 20/60	400 (4)	700 (7)	11
ABXL 22/40-1 200 ABXL 22/60-1 200	22/40 22/60	400 (4)	700 (7)	11
ABXL 24/40-1 200 ABXL 24/60-1 200	24/40 24/60	400 (4)	700 (7)	12
ABXL 26/40-1 200 ABXL 26/60-1 200	26/40 26/60	400 (4)	700 (7)	13
ABXL 28/40-1 200 ABXL 28/60-1 200	28/40 28/60	400 (4)	700 (7)	13
ABXL 30/40-1 200 ABXL 30/60-1 200	30/40 30/60	400 (4)	600 (6)	13
ABXL 32/40-1 200 ABXL 32/60-1 200	32/40 32/60	400 (4)	600 (6)	13

Packaging Unit: 1

Order Codes Andra2Balloon PTA-Catheter

Order Code	Innerballoon Ø/Length in mm	Outerballoon Ø/Length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
A2B 20/40-1 000 A2B 20/60-1 000	10/40 10/60	20/40 20/60	400 (4)	700 (7)	11
A2B 22/40-1 000 A2B 22/60-1 000	10/40 10/60	22/40 22/60	400 (4)	700 (7)	12
A2B 24/40-1 000 A2B 24/60-1 000	14/40 14/60	24/40 24/60	400 (4)	700 (7)	13
A2B 26/40-1 000 A2B 26/60-1 000	14/40 14/60	26/40 26/60	400 (4)	700 (7)	14
A2B 28/40-1 000 A2B 28/60-1 000	14/40 14/60	28/40 28/60	400 (4)	700 (7)	14
A2B 30/40-1 000 A2B 30/60-1 000	14/40 14/60	30/40 30/60	400 (4)	600 (6)	14
A2B 32/40-1 000 A2B 32/60-1 000	14/40 14/60	32/40 32/60	400 (4)	600 (6)	14

Packaging Unit: 1

Andramed GmbH

Schiesswiesenstr. 18
72766 Reutlingen, Germany
Tel.: +49 (0) 7127 811 57 15
Fax: +49 (0) 7127 811 57 11
info@andramed.com
www.andramed.com

Made in Germany
Rev. 03/2019

