


AndraStent® L/XL
Peripheral Bare Stent

Cobalt chromium bare metal stent
Diameters 6 – 25 mm
Lengths 13 – 57 mm

AndraStent[®] Peripheral Bare Stent L



Order Codes AndraStent[®] Peripheral Bare Stent L

Order Code	Length in mm	Expansion Stent Diameter in mm					Minimum recommended Diameter Introducer Sheath
		6	8	12	15	18	
AS13L	13	13	12	11	10	8	7 Fr*
AS17L	17	17	17	15	14	11	
AS21L	21	21	21	19	18	15	
AS26L	26	26	25	23	21	17	

* The diameter of the introducer sheath depends on the balloon (+ 2 Fr for the stent). Packaging Unit: 1

AndraStent® Peripheral Bare Stent XL



Order Codes AndraStent® Peripheral Bare Stent XL

Order Code	Length in mm	Expansion Stent Diameter in mm							Minimum recommended Diameter Introducer Sheath
		8	12	15	18	20	22	25	
AS13XL	13	13	12	12	11	10	9	8	8 Fr*
AS17XL	17	17	16	16	15	14	13	11	
AS21XL	21	21	20	20	19	18	16	14	
AS26XL	26	26	24	23	22	21	19	17	
AS30XL	30	30	29	28	27	25	23	23	
AS35XL	35	35	33	31	29	27	25	24	
AS39XL	39	39	36	35	33	31	29	27	
AS43XL	43	43	40	38	36	34	32	29	
AS48XL	48	48	44	42	40	37	35	33	
AS57XL	57	56	52	50	47	44	42	40	

* The diameter of the introducer sheath depends on the balloon (+ 2 Fr for the stent). Packaging Unit: 1

Andramed GmbH

Schiesswiesenstr. 18
72766 Reutlingen, Germany
Tel.: +49 (0) 7127 811 57 15
Fax: +49 (0) 7127 811 57 11
info@andramed.com
www.andramed.com

Made in Germany
Rev. 02/2019

